



LabQuip Limited  
The Black Barn,  
White Ox,  
Penrith,  
Cumbria.  
CA11 9NQ

Tel: +44 (0) 1768 895080  
E-mail: sales@lab-quip.co.uk  
Web-Site: www.lab-quip.co.uk  
VAT No: GB 698 4784 52  
Co Reg No: 033 99 679



## General Purpose Ovens 250°C to 300°C

### Features:

Temperature range 40°C to 250°C and 40°C to 300°C  
Aluminium coated mild steel or stainless steel interior  
Direct reading hydraulic controls  
Safety thermostat with indicator

### Optional extras:

Fan assisted circulation (30 to 350 litre only)  
Digital control and readout  
Programmable controllers  
Dial thermometer (30 to 350 litre only)  
Various timers (50 to 350 litre only)

### Design

The exterior is constructed from sheet steel finished in easy clean stove powder paint.  
The interior is fitted with fixed shelf runners with removable chrome plated steel wire grid shelves  
The door and chamber are fitted with silicone seals for excellent insulation against heat and is easy to clean  
The oven are heated by incoloy sheathed elements and are powered by 220/240v single phase supply

<b>GP Series</b>	<b>10</b>	<b>20</b>	<b>30</b>	<b>40</b>	<b>50</b>	<b>75</b>	<b>100</b>
Ext Dim HxWxD	415x375x350mm	540x375x382mm	570x415x422mm	670x415x422mm	615x590x422mm	623x590x567mm	735x590x567mm
Int Dim HxWxD	150x275x243mm	275x275x275mm	305x315x315mm	405x315x315mm	325x490x315mm	333x490x460mm	445x490x460mm
Volume (litres)	10	20	30	40	50	75	100
Shelf Positions	1	2	3	3	3	3	4
Number of shelves	1	2	2	2	2	2	3
Max load height	105mm	190mm	220mm	320mm	240mm	248mm	360mm
Shelf size WxD	240x208mm	240x240mm	280x280mm	280x280mm	455x270mm	455x415mm	455x415mm
Power rating (Watts)	350	500	750	1000	1000	1250	1250
Net weight	10kg	15kg	20kg	25kg	30kg	35kg	40kg
Fluctuation (hydraulic)	+/- 0.75°C						
Fluctuation (Digital)	+/- 0.5°C						

<b>GP Series</b>	<b>120</b>	<b>150</b>	<b>180</b>	<b>210</b>	<b>240</b>	<b>280</b>	<b>350</b>
Ext Dim HxWxD	735x590x659mm	845x590x659mm	957x590x659mm	957x590x749mm	1178x590x659mm	941x750x769mm	1104x750x769mm
Int Dim HxWxD	445x490x552mm	555x490x552mm	667x490x552mm	667x490x642mm	888x490x552mm	651x650x662mm	814x650x662mm
Volume (litres)	120	150	180	210	240	280	350
Shelf Positions	4	5	6	6	8	6	7
Number of shelves	3	3	4	4	4	4	4
Max load height	360mm	470mm	582mm	582mm	803mm	566mm	729mm
Shelf size WxD	455x507mm	455x507mm	455x507mm	455x597mm	455x507mm	615x617mm	615x617mm
Power rating (Watts)	1500	1500	1750	2000	2250	2250	2500
Net weight	45kg	50	55kg	65kg	73kg	78kg	95kg
Fluctuation (hydraulic)	+/- 0.75°C						
Fluctuation (Digital)	+/- 0.5°C						