

SieveMaster SQD



The SieveMaster SQD is a maintenance free shaker with a stainless steel body complete with aluminium top plate.

Capacity:

10 x 200mm dia sieves

8 x 300mm dia sieves.

Dimensions:

Body – 400mm x 400mm x 225mm

Top plate – 450mm x 450mm

Power:

240 volt single phase



The digital timer has a 99 minute countdown.

It is remote from the machine to allow fixing to an adjacent wall.

We can also supply this machine to take 400mm & 450mm dia sieves